

Mouse Anti-Human IgM (Fc) antibody [CH2] (HRP)

Cat. No. GTX20772-01

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human IgM.
Purification	Purified IgG
Conjugation	Horseradish peroxidase(HRP) horseradish peroxidase (HRP) RZ=3


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX20772-01 WB Image

WB analysis of human plasma using GTX20772-01 Human IgM (Fc) antibody [CH2] (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).