

Mouse Anti-Human IgM (Fc) antibody [CH2] (PerCP-Cy5.5)

Cat. No. GTX20772-11

| | |
|---------------------|------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Applications | FCM |
| Reactivity | Human |

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| FCM | Assay dependent |

Not tested in other applications.

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 15mM Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light. |
| Concentration | 0.1 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Purified human IgM. |
| Purification | Purified by size-exclusion chromatography |
| Conjugation | Peridinin-chlorophyll proteins-Cyanine5.5 (PerCP-Cy5.5) Wavelength |
| Note | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Cy and CyDye are trademarks of GE Healthcare. |



For full product information, images and publications, please visit our [website](#).