

FGFR1 antibody [M2F12]

Cat. No. GTX20829

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IHC-P, IHC-Fr, FCM, IP
Reactivity	Human, Mouse, Rat

References (4)

Package

125 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/ml
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent

Note : Recognizes 10 ng recombinant bacterial, baculo-viral, or native FGFr1.

Not tested in other applications.

Calculated MW 92 kDa. ([Note](#))**Product Note** This antibody reacts with the NH2-terminus of unique NH2-terminal Ig loop of FGFr1. Epitope is within the sequence between glu30 and ala74 of FGFr1α.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human ectodomain of FGFr1α expressed in E. coli beginning with pro23; antigen contained NH2-terminal gly-ser-pro-gly-ile and COOH- terminal glu-phe sequences.
Purification	Protein G purified

For full product information, images and publications, please visit our [website](#).

Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).