

HMW Kininogen antibody [C11C1]

Cat. No. GTX21004

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	Coagulation assay, WB, ELISA
Reactivity	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Coagulation assay	Assay dependent

Not tested in other applications.

Calculated MW 72 kDa. ([Note](#))

Product Note This antibody recognizes an epitope on the unique 46kDa light chain of human HMWK. The epitope maps to domain 5 of human HMWK. It can completely neutralize the coagulant activity of human HMWK.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human high molecular weight kininogen.
Purification	Protein G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX21004 WB Image

WB analysis of (A) normal human plasma, (B) HMWK-deficient plasma, and (C) purified HMWK proteins using GTX21004 HMW Kininogen antibody [C11C1].



For full product information, images and publications, please visit our [website](#).