

CD8 antibody [MEM-31] (PE-Cy7)

Cat. No. GTX21087-10

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM
Reactivity	Human

Package
100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	4 μ l reagent / 100 μ l of whole blood or 10^6 cells in a suspension

Not tested in other applications.

Product Note

This antibody recognizes a conformationally-dependent epitope of CD8. The antibody does not react with formaldehyde-fixed cells; negative in Western blotting application.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Crude thymus membrane fraction.
Purification	Purified by size-exclusion chromatography
Conjugation	Phycoerythrin-Cyanine7 (PE-Cy7) Wavelength
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Cy and CyDye are trademarks of GE Healthcare.



For full product information, images and publications, please visit our [website](#).