

Multiple Species Tissue Array II - Normal Tissues

Cat. No. GTX21333

Application IHC-P, ISH

Package
5 slide

PRODUCT

Summary

GeneTex's Multiple Species Paraffin Tissue Array is designed for testing the species reactivity of antibodies. This array can bring very important information of antibodies in immunohistochemistry studies. The same type of tissues from 15 different species, including human, mouse, rat, Rhesus monkey, Cynomolgus monkey, rabbit, guinea pig, goat, sheep, hamster, chicken, bovine, horse, dog, and porcine, are arrayed on positively charged glass slides. The Multiple Species paraffin tissue array can also enable you to study expression of genes or proteins of interest among different animals.

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
ISH	Assay dependent

Not tested in other applications.

PROPERTIES

Storage Cool dry at room temperature. 4°C for long term storage.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).