

Cortactin antibody, N-term

Cat. No. GTX21374

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.01-0.03µg/ml

Not tested in other applications.

Calculated MW 62 kDa. ([Note](#))**Product Note** This antibody is expected to recognize both reported isoforms (NP_005222.2; NP_612632.1).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence WKASAGHAVSIA-C, from the N Terminus of the protein sequence according to NP_005222.2; NP_612632.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES

GTX21374 WB Image

WB analysis of human brain (cerebellum) lysate using GTX21374 Cortactin antibody, N-term.

Dilution : 0.1 μ g/ml


Loading : 35 μ g protein in RIPA buffer


GTX21374 WB Image

WB analysis of HEK293 cell lysate using GTX21374 Cortactin antibody, N-term.

Dilution : 0.03 μ g/ml

Loading : 35 μ g protein in RIPA buffer


GTX21374 WB Image

WB analysis of mouse brain lysate using GTX21374 Cortactin antibody, N-term.

Dilution : 0.01 μ g/ml

Loading : 35 μ g protein in RIPA buffer



For full product information, images and publications, please visit our [website](https://www.genetex.com).