

Human Small Intestine: Ileum (Diabetes) tissue slides

Cat. No. GTX21798

Applications	IHC-P, ISH
Species	Human

Package
5 slide

PRODUCT

Summary

Single human diabetes diseased tissue with 5 µm thickness is mounted on each positively charged glass slide. These tissue sections are designed for rapidly identifying cellular localization of genes or proteins in human diabetes diseased tissues. More organs may be available upon request.

Tissue type

Diabetes

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution**Recommended dilution**

IHC-P	Assay dependent
ISH	Assay dependent

Not tested in other applications.

Properties

Storage

Cool dry at room temperature. 4°C for long term storage.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).