

## Human Skeletal Muscle (Diabetes) tissue slides

Cat. No. GTX21805

Applications IHC-P, ISH

Species Human

References ( 1 )

Package

5 slide

## PRODUCT

## Summary

Single human diabetes diseased tissue with 5  $\mu$ m thickness is mounted on each positively charged glass slide. These tissue sections are designed for rapidly identifying cellular localization of genes or proteins in human diabetes diseased tissues. More organs may be available upon request.

## Tissue type

Diabetes

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

## Suggested dilution

## Recommended dilution

IHC-P

Assay dependent

ISH

Assay dependent

Not tested in other applications.

## Properties

## Storage


Cool dry at room temperature. 4°C for long term storage.

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



GTX21805 IHC-P Image



For full product information, images and publications, please visit our [website](#).