

Bacillus anthracis spore antigen antibody [SA27]

Cat. No. GTX21994

| | |
|--------------|--------------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG2a |
| Applications | WB, ELISA |
| Reactivity | Bacillus anthracis |

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB | Assay dependent |
| ELISA | Assay dependent |

Note : This antibody could be used as detection antibody.

Not tested in other applications.

Product Note No cross-reactivity found in ELISA with B. anthracis vegetative cells, B. globigii, B. subtilis and B. cereus spores.

Properties

| | |
|---------------|-----------------------------------------------------------------------------------------|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 0.09% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
| Immunogen | B. anthracis spore extract |
| Purification | Protein G purified |
| Conjugation | Unconjugated |

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).