

Human Eye: Retina (Normal) tissue slides

Cat. No. GTX22016

Application IHC-P, ISH

Species Human

★★★★★ Review (1)

Package
5 slide

PRODUCT

Summary

A single human adult normal tissue with 5 µm thickness is mounted on each positively charged glass slide. These tissue sections are designed for rapidly identifying cellular localization of genes or proteins in human normal tissues. More organs maybe available upon request.

Tissue type Normal

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
ISH	Assay dependent

Not tested in other applications.

PROPERTIES

Storage Cool dry at room temperature. 4°C for long term storage.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).