

GPBB antibody [M8030622]

Cat. No. GTX22075

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

ELISA: Use at an assay dependent dilution. Functions as a matched pair (conjugate) with GTX22069 as the solid phase. Optimal dilutions/concentrations should be determined by the end user.

Product Note Affinity Constant: 1.0×10^{10} L/M

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	10.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human Glycogen phosphorylase isoenzyme BB
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).