

## Cyclin E antibody

## Cat. No. GTX22094

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, ELISA
Reactivity	Human

Package 100 μl

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

^Optimal dilutions/con	ncentrations should be determined by the researcher.	
Suggested dilution	Recommended dilution	
WB	1:2000-1:10000	
IP	1:100	
ELISA	1:10000	
Not tested in other app	plications.	
Calculated MW	47 kDa. ( <u>Note</u> )	
Product Note	No reaction was observed against other related cyclins.	
Properties		
_		

Properties	
Form	Liquid
Buffer	Filtered-steriled serum
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	A KLH conjugated peptide corresponding to a region near the carboxy terminus of human Cyclin E.
Purification	Unpurified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 1