

Cyclin T1 antibody

Cat. No. GTX22098

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
IHC-P	1:200-1:1000
IP	1:100
ELISA	1:1000-1:5000

Not tested in other applications.

Calculated MW 81 kDa. ([Note](#))**Product Note** No reaction was observed against other related cyclins.

Properties

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	75 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Cyclin T1 peptide corresponding to an internal region of the human protein conjugated to KLH.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX22098 IHC-P Image

IHC-P analysis of human skin tissue using GTX22098 Cyclin T1 antibody.

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).