

## Human Brain: Medulla oblongata (Parkinson's Disease) tissue slides

Cat. No. GTX22155

**Application** IHC-P, ISH**Species** Human**Package**

5 slide

## PRODUCT

**Summary**

Single human neurological diseased tissue with 5  $\mu$ m thickness is mounted on each positively charged glass slide. These tissue sections are designed for rapidly identifying cellular localization of genes or proteins in human Alzheimer, Parkinson's Disease, Dementia, Multiple Sclerosis, Progressive Supranuclear Palsy (PSP), and Depression diseased tissues. More organs maybe available upon request.

**Tissue type** Neurological Disease

## APPLICATION

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

**Suggested dilution****Recommended dilution**

IHC-P

Assay dependent

ISH

Assay dependent

Not tested in other applications.

## PROPERTIES

**Storage**

Cool dry at room temperature. 4°C for long term storage.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).