

EasyBlot anti Mouse IgG (HRP)

Cat. No. GTX221667-01

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Mouse

References (39)

★★★★★ Review (1)

Package

250 µl, 50 µl

PRODUCT

Summary

EasyBlot anti-Mouse IgG is an HRP-conjugated second step reagent that specifically reacts with the native, non-reduced form of mouse IgG and does not bind the reduced, denatured form. Using Easyblot decreases the interference caused by the heavy (~50 kDa) and light chains (~25 kDa) of the IgG used for immunoprecipitation (IP). EasyBlot is simple to use and results in Western Blot detection of only the target band in your IP experiments. Exercise caution when using Protein A/G for IP with EasyBlot. Protein A/G contaminants bind to mouse IgG with high affinity and may cause contaminating bands with a molecular weight similar to IgG heavy chain or light chain.

To minimize background from Protein A/G contamination, we strongly recommend the use of EasyBlot Kit anti-Mouse IgG (Optimized for Prot A/G) (GTXGTX225857-01) or EasyBlokter (GTX425858) to minimize protein A/G contamination or you can use anti-Mouse IgG Agarose beads (GTX27018) instead of Protein A/G.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution Recommended dilution

WB	1:1000
----	--------

Note : Use 2-5 µg of the IP antibody; do not exceed 5 µg. Completely reduce all IP samples by 100mM DTT (fresh prepared). For further detail, please refer to the protocol.

Not tested in other applications.

Properties

Form	Liquid
Buffer	0.05M Tris, 0.15M NaCl, 1%BSA, 40% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.03 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified by antigen-affinity chromatography.
Conjugation	Horseradish peroxidase(HRP)



For full product information, images and publications, please visit our [website](#).


Date 2025 / 05 / 27 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx221667-01 WB Image

Western Blot analysis of native mouse IgG (left) or denatured mouse IgG (right) using EasyBlot anti Mouse IgG (HRP) antibody.



GTx221667-01 WB Image

EasyBlot anti-Mouse IgG (HRP) (GTx221667-01) was used as a secondary reagent for immunoprecipitation analysis. 293T lysates were immunoprecipitated using 5 ug of normal mouse IgG or ELP3 antibody [GT11811] (GTx631395)



For full product information, images and publications, please visit our [website](#).