

ELMO2 antibody

Cat. No. GTX22240

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

Package

50 µg

Applications

Application Note

Western Blot - Approx 80kDa band seen in Human Brain lysate. Recommended for use at 1.0-3.0µg/ml.

Calculated MW	83 kDa. (Note)
Product Note	This antibody is expected to recognise both human isoforms.


Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence PKEPSSYDFVYHYG, from C Terminus of the protein sequence according to NP_071369; NP_573403
Purification	IgG fraction Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX22240 WB Image

Western blot analysis of ELMO in human brain lysate using GTX22240 at 2 ug/ ml.



GTX22240 WB Image



For full product information, images and publications, please visit our [website](https://www.genetex.com).