

# 12L HETE antibody

**Cat. No. GTX22275**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	ELISA, RIA
<b>Reactivity</b>	Species independent

**Package**  
100 µl

## Applications

### Application Note

Working dilutions are best determined by the end user.

**Product Note** This antibody cross reacts with Methyl ester-12-HETE (20%), 8-HETE (8.7%) and 12(L)-HHT (8.6%).

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	unknown
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Immunogen</b>	Pure 12(L)-HETE conjugated to human serum albumin
<b>Purification</b>	Whole antiserum
<b>Conjugation</b>	Unconjugated
<b>Note</b>	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.  Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).