

12L HETE antibody

Cat. No. GTX22275

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	ELISA, RIA
Reactivity	Species independent

Package
100 µl

APPLICATION

Application Note

Working dilutions are best determined by the end user.

Product Note This antibody cross reacts with Methyl ester-12-HETE (20%), 8-HETE (8.7%) and 12(L)-HHT (8.6%).

PROPERTIES

Form	Liquid
Buffer	unknown
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Pure 12(L)-HETE conjugated to human serum albumin
Purification	Whole antiserum
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).