

15 deoxy Delta 12,14 Prostaglandin J2 antibody

Cat. No. GTX22313

| | |
|---------------------|---------------------|
| Host | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Applications | ELISA, RIA |
| Reactivity | Species independent |

Package

50 µl

Applications

Application Note

Working dilutions are best determined by the end user.

Product Note

This antibody detects Prostaglandin J2 (25.3 %), Prostaglandin D2 (13.13 %), D12-Prostaglandin J2 (8.89 %) and 9a,11b-Prostaglandin F2 alpha (6.93 %).

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 0.1% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | 15-deoxy-D12,14-Prostaglandin J2 conjugated to a protein carrier molecule |
| Conjugation | Unconjugated |

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).