

PIP2 antibody

Cat. No. GTX22335

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ELISA, RIA
Reactivity	Species independent

Package
50 µg

Applications

Application Note

Working dilutions are best determined by the end user.

Product Note This antibody cross reacts with Phosphatidic acid (<5.0%).

Properties

Form	Liquid
Buffer	PBS, 10% newborn calf Serum
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Phosphatidylinositol-4,5-bisphosphate (PIP2) from bovine spinal cord
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).