

Her2 / ErbB2 peptide

Cat. No. GTX22429

Application	IP
Species	Human

Package
500 µl

APPLICATION

Application Note

Blocking peptide for rabbit polyclonal antibody to Her2 / ErbB2(GTX22428).

Calculated MW 138 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS, BSA
Preservative	0.1% Sodium azide
Storage	Store this product at 4°C, DO NOT FREEZE. The product is stable for one year from the date of shipment.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	Peptide mapping to the carboxy terminus of human Neu (HER2 or ErbB-2).
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).