

# APAF1 peptide

**Cat. No. GTX22446**

<b>Applications</b>	IP
<b>Species</b>	Human

**Package**  
500 µl

## Applications

### Application Note

Blocking peptide for rabbit polyclonal antibody to APAF1(GTX22445).

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, BSA
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store this product at 4°C, DO NOT FREEZE. The product is stable for one year from the date of shipment.
<b>Concentration</b>	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Region/Sequence</b>	Peptide mapping to the amino terminus of human Apaf-1.
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).