

p107 peptide

Cat. No. GTX22452

| | |
|--------------------|-------|
| Application | IP |
| Species | Human |

Package
500 µl

APPLICATION

Application Note

Blocking peptide for rabbit polyclonal antibody to p107(GTX22451).

Calculated MW 121 kDa. ([Note](#))

PROPERTIES

| | |
|------------------------|---|
| Form | Liquid |
| Buffer | PBS, BSA |
| Preservative | 0.1% Sodium azide |
| Storage | Store this product at 4°C, DO NOT FREEZE. The product is stable for one year from the date of shipment. |
| Concentration | 0.2 mg/ml (Please refer to the vial label for the specific concentration.) |
| Region/Sequence | Peptide mapping to the carboxy terminal domain of human p107. |
| Conjugation | Unconjugated |

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).