

CD38 antibody [HIT2] (PE-Cy5)

Cat. No. GTX22577-09

| | |
|---------------------|------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Applications | FCM |
| Reactivity | Human |

Package
100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|---|
| FCM | 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension |

Not tested in other applications.

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 15mM Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
| Immunogen | Human thymocytes in foetus |
| Purification | Purified by size-exclusion chromatography |
| Conjugation | Phycoerythrin-Cyanine5 (PE-Cy5) Wavelength |
| Note | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Cy and CyDye are trademarks of GE Healthcare. |



For full product information, images and publications, please visit our [website](#).