

## Hepatitis C virus Core Antigen antibody [1E5] (FITC)

**Cat. No. GTX22584**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	IP, ELISA
<b>Reactivity</b>	Hepatitis C virus

**Package**  
50 µg

## Applications

**Application Note**

Optimal working dilutions are best determined by the investigator.

**Product Note** This antibody is specific for Hepatitis C Core Antigen.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.01% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant fragment, corresponding to amino acids 1-80 of Hepatitis C Core Antigen.
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).