

## EasyBlot anti-Mouse IgG Kit (Optimized for Prot A/G)

Cat. No. GTX225857-01

Applications	WB
Reactivity	Mouse

References ( 4 )  
Package  
250 µl, 50 µl

## PRODUCT

Content	EasyBlot Kit Package	Content	Package	Storage Temp
	GTX225857-01 (250µl)	GTX221667-01 EasyBlot anti-Mouse IgG (HRP)	250µl	-20 °C
		GTX425858 EasyBlocker	60g	4 °C
	GTX225857-01 (50µl)	GTX221667-01 EasyBlot anti-Mouse IgG (HRP)	50µl	-20 °C
		GTX425858 EasyBlocker	12g	4 °C

## Summary

EasyBlot anti-Mouse IgG is an HRP-conjugated second step reagent that specifically reacts with the native, non-reduced form of mouse IgG and does not bind to the reduced, denatured form. Using Easyblot decreases the interference caused by the heavy (~50 kDa) and light chains (~25 kDa) of the IgG used for immunoprecipitation (IP).

The protein A/G beads used in the IP procedure sometimes bind to IgG and produce contaminating signals at approximately 50 kDa and 25 kDa. This kit also includes an EasyBlocker (GTX425858) that can minimize the background caused by Protein A/G contamination. GeneTex's EasyBlot kit is simple to use and detects only the target band in your IP experiments, resulting in cleaner western blots.

## Applications

## Application Note

Suggested starting dilutions for WB is 1:1000. Use 2-5 µg of the IP antibody; do not exceed 5 µg. Completely reduce all IP samples by 100mM DTT (fresh prepared). For further detail, please refer to the protocol. Not yet tested in other applications. Optimal working dilutions should be determined experimentally by the end user.

## Properties

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.03 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified by antigen-affinity chromatography.
Conjugation	Horseradish peroxidase(HRP)



For full product information, images and publications, please visit our [website](#).


Date 2025 / 12 / 19 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX225857-01 WB Image

EasyBlot anti-Mouse IgG Kit (Optimized for Prot A/G) (GTX225857-01) was used as a secondary reagent for immunoprecipitation analysis.

HepG2 lysates were immunoprecipitated using 5 ug of normal rabbit IgG or ELP3 antibody [GT11811] (GTX631395).



For full product information, images and publications, please visit our [website](https://www.genetex.com).