

PAP pen

Cat. No. GTX22601

Applications IHC-P, IHC-Fr, ISH**Package**
1 unit

PRODUCT

Summary PAP hydrophobic barrier pen for immunohistochemistry

Applications

Application Note

The PAP pen is a special marking pen that provides a thin film-like green-tinged hydrophobic barrier when a circle is drawn around a specimen on a slide. This water repellent barrier keeps staining reagents localized on the tissue sections, prevents mixing of reagents when differentially staining two sections on the same slide, and allows use of less reagents per section. Insoluble in alcohol and acetone, PAP Pen markings can be optionally removed with xylene after the staining procedure is complete. PAP Pen is good for approximately 500 uses.

Properties

Storage Store at room temperature.**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX22601 Image



For full product information, images and publications, please visit our [website](#).

Date 2025 / 08 / 02 Page 1 of 1