

## BrdU antibody [BrdU007] (ready-to-use) (Biotin)

**Cat. No. GTX22704**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	ICC/IF, IHC-P, IHC
<b>Reactivity</b>	Species independent

**Package**  
6 ml

## Applications

**Application Note**

This antibody is ready to use and has been optimised for use in IHC, so that it needs no dilution.

**Product Note** This antibody recognizes bromodeoxuridine (BrdU).

## Properties


<b>Form</b>	Liquid
<b>Buffer</b>	Tissue culture supernatant diluted with PBS, protein carrier
<b>Preservative</b>	Proprietary preservative , not Sodium azide or Thimerosal
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Conjugation</b>	Biotin

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES

**GTX22704 IHC Image**

Bromodeoxyuridine positive mouse intestine stained with anti-BrdU antibody (GTX22704).



For full product information, images and publications, please visit our [website](#).