

BrdU antibody [BrdU007] (ready-to-use) (Biotin)

Cat. No. GTX22704

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, IHC-P, IHC
Reactivity	Species independent

Package
6 ml

Applications

Application Note

This antibody is ready to use and has been optimised for use in IHC, so that it needs no dilution.

Product Note This antibody recognizes bromodeoxuridine (BrdU).

Properties


Form	Liquid
Buffer	Tissue culture supernatant diluted with PBS, protein carrier
Preservative	Proprietary preservative , not Sodium azide or Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Conjugation	Biotin

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX22704 IHC Image**

Bromodeoxyuridine positive mouse intestine stained with anti-BrdU antibody (GTX22704).



For full product information, images and publications, please visit our [website](#).