

EasyBlot anti-Rat IgG Kit (Optimized for Prot A/G)

Cat. No. GTX228909-01

| | |
|---------------------|-----|
| Applications | WB |
| Reactivity | Rat |

Package
250 µl, 50 µl

PRODUCT

| | EasyBlot Kit Package | Content | Package | Storage Temp |
|----------------|----------------------|--|---------|--------------|
| Content | GTX228909-01 (250µl) | GTX628474-01 EasyBlot anti Rat IgG (HRP) | 250µl | -20 °C |
| | | GTX425858 EasyBlocker | 60g | 4 °C |
| | GTX228909-01 (50µl) | GTX628474-01 EasyBlot anti Rat IgG (HRP) | 50µl | -20 °C |
| | | GTX425858 EasyBlocker | 12g | 4 °C |

EasyBlot anti-Rat IgG is an HRP-conjugated second step reagent that specifically reacts with the native, non-reduced form of rat IgG and does not bind to the reduced, denatured form. Using Easyblot decreases the interference caused by the heavy (~50 kDa) and light chains (~25 kDa) of the IgG used for immunoprecipitation (IP).

Summary

The protein A/G beads used in the IP procedure sometimes bind to IgG and produce contaminating signals at approximately 50 kDa and 25 kDa. This kit also includes an EasyBlocker (GTX425858) that can minimize the background caused by Protein A/G contamination. GeneTex's EasyBlot kit is simple to use and detects only the target band in your IP experiments, resulting in cleaner western blots.

Applications

Application Note

Suggested starting dilutions for WB is 1:1000. Use 2-5 µg of the IP antibody; do not exceed 5 µg. Completely reduce all IP samples by 100mM DTT (fresh prepared). For further detail, please refer to the protocol. Not yet tested in other applications. Optimal working dilutions should be determined experimentally by the end user.

Properties

| | |
|----------------------|--|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 0.37 mg/ml (Please refer to the vial label for the specific concentration.) |
| Purification | Affinity purified by Protein G. |
| Conjugation | Horseradish peroxidase(HRP) |



For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).