

CDK8 antibody

Cat. No. GTX22955

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human, Mouse, Rat

Package
500 µl

Applications

Application Note

Immunofluorescence: assay dependent. Immunoprecipitation (with protein A, native): use at a concentration of 10µg / mg of lysate. Western Blotting: use at a concentration of 1 - 2µg / ml for 2hrs at RT. Optimal dilutions / concentrations should be determined by the end user. By Western blot, this antibody detects a band of 54 kDa, which corresponds to the predicted molecular weight of Cdk 8.

Calculated MW 53 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide: C-SQQPPQYSHQTHRY, corresponding to amino acids 451-464 of Human CDK8.
Purification	Protein A affinity purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 1