

## STAT4 antibody [ST4.01 (NB.34)]

Cat. No. GTX23285

Host	Armenian Hamster
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IP, Gel supershift assays
Reactivity	Human, Mouse

Package  
500 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent
Gel supershift assays	Assay dependent

Not tested in other applications.

Calculated MW 86 kDa. ( [Note](#) )

**Product Note** Please note that this MAb is of Armenian Hamster origin and use the appropriate secondary Ab (from Jackson ImmunoResearch) for all applications.

## Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length protein (Mouse).
Purification	Protein G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).