

LMP2 antibody

Cat. No. GTX23328

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-3 µg/ml
Not tested in other applications.	

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mixture of synthetic peptides corresponding to residues C R(207) V I L G N/D E L P K F Y D E(220) of human and mouse proteasome 20S LMP2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX23328 WB Image

WB analysis of 50 µg mouse embryo fibroblast cell extract using GTX23328 LMP2 antibody.


GTX23328 WB Image

WB analysis of whole cell extracts (30 µg lysate) of HeLa (Lane 1), HeLa treated with IFN gamma (100ng/ml IFN gamma for 72h) (Lane 2), THP-1 (Lane 3), THP-1 treated with PMA (50 ng PMA for 48h) and LPS (200 ng/ml LPS for 24 h) (Lane 4), U-937 (Lane 5), Raji (Lane 6) and Ramos (Lane 7) using GTX23328 LMP2 antibody.

Dilution : 2 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).