

MRP1 antibody [QCRL3]

Cat. No. GTX23370

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	IHC
Reactivity	Human

Package
250 µl

Applications

Application Note

IHC: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

Product Note

This antibody detects a 190 kDa transmembrane phosphoglycoprotein overexpressed in various humation dependent multidrug resistant tumor cell lines. It does not cross react with human MDR1 and MDR3, or mouse MRP gene products.

Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Membranes prepared from the human small cell lung cancer cell line H69AR which highly overexpresses MRP
Purification	IgG fraction
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).