

# MRP2 antibody [M2III-6]

**Cat. No. GTX23373**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	WB, ICC/IF, IHC-P, IHC-Fr
<b>Reactivity</b>	Human, Rat

**Package**  
250 µl

## Applications

### Application Note

Optimal dilutions / concentrations should be determined by the end user.

<b>Calculated MW</b>	174 kDa. ( <a href="#">Note</a> )
<b>Product Note</b>	This antibody detects MRP 2. It does not cross-react with human MDR 1, MRP 1, MRP 3 or MRP 5 gene products.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 1% BSA
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage 1-2 weeks, store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Bacterial fusion protein of human MRP 2 containing the carboxy-terminal region (202-)
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).