

MRP2 antibody [M2III-6]

Cat. No. GTX23373

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-P, IHC-Fr
Reactivity	Human, Rat

Package
250 µl

Applications

Application Note

Optimal dilutions / concentrations should be determined by the end user.

Calculated MW	174 kDa. (Note)
Product Note	This antibody detects MRP 2. It does not cross-react with human MDR 1, MRP 1, MRP 3 or MRP 5 gene products.

Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage 1-2 weeks, store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Bacterial fusion protein of human MRP 2 containing the carboxy-terminal region (202-)
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).