

ABCG2 antibody [BXP-34]

Cat. No. GTX23379

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	ICC/IF, IHC-Fr
Reactivity	Human

Package
250 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-Fr	Assay dependent

Not tested in other applications.

Product Note This antibody does not cross react with the human MDR 1, MRP 1, MRP 2, or MRP 5 gene products.

Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Immunogen	BCRP over-expressing cell line MCF 7 MR.
Purification	Tissue Culture Supernatant From Tissue Culture Supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).