

Cannabinoid receptor 2 antibody

Cat. No. GTX23561

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM, IHC
Reactivity	Human, Mouse, Rat

References (4)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:50-1:500
ICC/IF	1:500
IHC-P	1:10-1:100
FCM	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 39 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS, 0.1% BSA, 50% Glycerol
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A fusion protein containing the first 32 amino acid residues from rat CB2.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTx23561 ICC/IF Image

ICC/IF analysis of AtT20 cells transfected with the rat CB2 using GTx23561 Cannabinoid receptor 2 antibody.


GTx23561 WB Image

WB analysis of 25 ug of HT29 (lane 1), C6 (lane 2) and rat colon (lane 3) cell lysates using GTx23561 Cannabinoid receptor 2 antibody.

Dilution : 1:200


GTx23561 IHC-P Image


IHC-P analysis of human tonsil tissue using GTx23561 Cannabinoid receptor 2 antibody.

Right : Primary antibody

Left : Negative control without primary antibody

Antigen retrieval : 10mM sodium citrate (pH 6.0) microwaved for 8-15 min

Dilution : 1:20


GTx23561 IHC-P Image

IHC-P analysis of human skin tissue using GTx23561 Cannabinoid receptor 2 antibody.

Right : Primary antibody

Left : Negative control without primary antibody

Antigen retrieval : 10mM sodium citrate (pH 6.0) microwaved for 8-15 min

Dilution : 1:20



For full product information, images and publications, please visit our [website](https://www.genetex.com).