

AMPK alpha 1 antibody

Cat. No. GTX23759

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, IHC
Reactivity	Human, Mouse, Rat, Bovine

References (1)

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2µg/ml
IP	Assay dependent
IHC	0.2- 2µg/ml

Not tested in other applications.

Calculated MW 64 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Tris-Citrate/Phosphate
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide representing a portion of the protein encoded within exon 8.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX23759 WB Image

WB analysis of (1) bovine aortic endothelial cell, (2) rat aortic smooth muscle cell, (3) HepG2 cell, and (4) human aortic endothelial cell lysates using GTX23759 AMPK alpha 1 antibody.

Dilution : 2µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).