

HEK293 whole cell lysate

Cat. No. GTX23872

Species Human

Package
100 µg

PRODUCT

Summary

HEK293 cells from ATCC (Cat# CRL-1573) were grown to 85% confluency, harvested, homogenized in modified RIPA buffer and cleared by centrifugation. Modified RIPA buffer. Protein concentration of the cell lysate was determined, and the lysate was boiled for 5 min in SDS sample buffer.

Properties

Form Liquid

Buffer 1x SDS sample buffer, 5% β-mercaptoethanol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 2 mg/ml (Please refer to the vial label for the specific concentration.)

Purification Whole Cell Lysate

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).