

SK-N-SH whole cell lysate

Cat. No. GTX23956

Species	Human
---------	-------

Package	100 µg
---------	--------

PRODUCT

Summary

SK-N-SH lysate was prepared by homogenization in modified RIPA buffer. Cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. The lysate was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Properties

Form	Liquid
------	--------

Buffer	1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	---

Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).

Date 2026 / 01 / 15 Page 1 of 1