

Blood Group Lewis a antibody [7LE]

Cat. No. GTX23967

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mucin isolated by size-exclusion chromatography from an ovarian cyst fluid (pure endocervical type according to the Fenoglio's classification).
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX23967 IHC-P Image

IHC-P analysis of colorectal adenocarcinoma tissue using GTX23967 Blood Group Lewis a antibody [7LE].



GTX23967 IHC-P Image

IHC-P analysis of human pancreas tissue using GTX23967 Blood Group Lewis a antibody [7LE].



GTX23967 IHC-P Image

IHC-P analysis of human colon adenocarcinoma tissue using GTX23967 Blood Group Lewis a antibody [7LE].



For full product information, images and publications, please visit our [website](#).