

Rat Pancreas (Normal) tissue lysate

Cat. No. GTX24036
Species Rat

Package
100 µg

PRODUCT
Summary

Pancreas lysate was prepared by homogenization in lysis buffer. Tissue and cell debris was removed by centrifugation. The product was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Tissue type

Normal

Properties
Form

Liquid

Buffer

1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol

Preservative

No preservatives

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration

2 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).