

## Ezrin antibody [3C12]

Cat. No. GTX24069

|              |              |
|--------------|--------------|
| Host         | Mouse        |
| Clonality    | Monoclonal   |
| Isotype      | IgG1         |
| Applications | IHC-P, IHC   |
| Reactivity   | Human, Mouse |

References ( 1 )

Package

500 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-P              | 1:50-1:100           |
| IHC                | Assay dependent      |

**Note : Recommendations for antigen retrieval : Citrate Buffer pH 6.0**

Not tested in other applications.

## Properties

|               |                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                     |
| Buffer        | PBS                                                                                                        |
| Preservative  | 0.09% Sodium azide                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.                    |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                           |
| Immunogen     | BALB/C mice were immunized with a recombinant fragment of human ezrin corresponding to amino acid 362-585. |
| Purification  | Purified IgG                                                                                               |
| Conjugation   | Unconjugated                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX24069 IHC-P Image

IHC-P analysis of human lung tissue using GTX24069 Ezrin antibody [3C12].



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 2 of 2