

## Rabbit Anti-Human kappa light chain antibody (FITC)

Cat. No. GTX24208

|              |               |
|--------------|---------------|
| Host         | Rabbit        |
| Clonality    | Polyclonal    |
| Isotype      | IgG           |
| Applications | ICC/IF, IHC-P |
| Reactivity   | Human         |

Package  
500 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ICC/IF             | Assay dependent      |
| IHC-P              | 1:20-1:40            |

**Note : Antigen retrieval : Proteinase K 30 min at 37°C**

Not tested in other applications.

## Properties

|               |                                                                                                             |
|---------------|-------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                      |
| Buffer        | PBS                                                                                                         |
| Preservative  | 0.09% Sodium azide                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                            |
| Immunogen     | Pool of human kappa light chain isolated from the urine of patients with Bence Jones Proteinuria.           |
| Purification  | Purified IgG                                                                                                |
| Conjugation   | Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>                                                |

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).