

## Parietal lobe (human): normal tissue slides

**Cat. No. GTX24313**

|                     |            |
|---------------------|------------|
| <b>Applications</b> | IHC-P, ISH |
| <b>Species</b>      | Human      |

**Package**  
5 slide

## PRODUCT

**Summary** This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.

**Tissue type** Normal

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-P              | Assay dependent      |
| ISH                | Assay dependent      |

Not tested in other applications.

## Properties

**Storage** Cool dry at room temperature. 4°C for long term storage.

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note** Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).