

Human Skeletal muscle (Normal) tissue slides

Cat. No. GTX24366

Application	IHC-P, ISH
Species	Human

Reference (1)
Package
5 slide

PRODUCT

Summary

A single human adult normal tissue with 5 μ m thickness is mounted on each positively charged glass slide. These tissue sections are designed for rapidly identifying cellular localization of genes or proteins in human normal tissues. More organs may be available upon request.

Tissue type Normal

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

IHC-P	Assay dependent
-------	-----------------

ISH	Assay dependent
-----	-----------------

Not tested in other applications.

PROPERTIES


Storage Cool dry at room temperature. 4°C for long term storage.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX24366 IHC-P Image



For full product information, images and publications, please visit our [website](#).