

Human Kidney (Liver Cirrhosis) tissue slides

Cat. No. GTX24593

Applications IHC-P, ISH**Species** Human

Package

5 slide

PRODUCT

Summary

Single human tissue from liver cirrhosis patients is mounted on each positively charged glass slide. These 5 µm thickness tissue sections are designed for rapidly identifying cellular localization of genes or proteins in human Liver Cirrhosis diseased tissues. More organs may be available upon request.

Tissue type

Liver Cirrhosis

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution**Recommended dilution**

IHC-P

Assay dependent

ISH

Assay dependent

Not tested in other applications.

Properties

Storage

Cool dry at room temperature. 4°C for long term storage.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).