

Nop2p antibody [22G1]

Cat. No. GTX24638

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	WB, ICC/IF
Reactivity	Yeast

Package
250 µl

Applications

Application Note

WB: Use at a dilution of 1/10000 (cell lysates) - 1/25000 (nuclear fractions). For other (non-ECL) detection methods, use at a dilution of 1/1000 - 1/5000. Predicted molecular weight: 100 kDa. IF: Use at a dilution of 1/500 - 1/5000 (yeast cells). Optimal dilutions/concentrations should be determined by the end user.

Product Note

Specificity was determined by western blot analysis of extracts from a yeast strain expressing Nop2p under GAL-UAS control and from strains containing nop2 ts alleles that express Nop2p isoforms of altered molecular weight.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	10mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	This was raised against a yeast nucleolar prep, and screened by immunocytochemistry and western blotting.
Purification	Tissue culture supernatant Sterile filtered.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).